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



3	定量分析	將偵測到的元素訊號轉換成元素含量。	<table border="1"> <thead> <tr> <th>Name</th> <th>Peak BE</th> <th>At. %</th> </tr> </thead> <tbody> <tr> <td>Si2p</td> <td>99.22</td> <td>68.78</td> </tr> <tr> <td>O1s</td> <td>531.98</td> <td>19.12</td> </tr> <tr> <td>N1s</td> <td>396.86</td> <td>12.1</td> </tr> </tbody> </table>	Name	Peak BE	At. %	Si2p	99.22	68.78	O1s	531.98	19.12	N1s	396.86	12.1
Name	Peak BE	At. %													
Si2p	99.22	68.78													
O1s	531.98	19.12													
N1s	396.86	12.1													
4	XPS Depth Profile (縱深成份分析)	藉由控制離子束濺蝕樣品的表面，分析來自不同深度的電子訊號，而得到元素成份的縱深分佈圖。	 <p>The graph shows the atomic percent of four elements as a function of etch time. The Si 2p (Element) signal (red) increases from 0% at 0s to approximately 100% at 750s. The Si 2p (NISI) signal (purple) peaks at about 60% around 180s. The Ni 2p (NISI) signal (blue) peaks at about 30% around 340s. The O1s signal (green) peaks at about 60% around 30s and then decreases to near 0% by 180s.</p>												
5	XPS Mapping (成像分析)	偵測樣品上欲分析區域的元素訊號，得到分析區域表面元素之位置分佈影像圖。	 <p>The figure displays four XPS mapping images of a sample surface. The top-left image shows Hf4f signal (red) with circular features. The top-right image shows Ta4f Metal signal (blue) with circular features. The bottom-left image shows Ta4f Oxide signal (green) with circular features. The bottom-right image shows N1s signal (cyan) with a more uniform distribution. A 1000 μm scale bar is present in the bottom-right image.</p>												
6	XPS Line Scan (線掃描分析)	偵測樣品上欲分析之直線距離的元素訊號，得到線上表面元素的分佈圖。	 <p>The line scan graph plots intensity against distance in micrometers. The x-axis ranges from 0 to 200 μm. The y-axis is labeled 'Intensity'. The graph shows several distinct peaks at various distances, indicating the presence of different elements along the scan line.</p>												
7	Angle-resolved Analysis (ARXPS) (角度解析分析)	可同時檢測不同角度之薄膜樣品 (厚度<7nm)表面成份鍵結分析。	 <p>The ARXPS plot shows binding energy in eV on the x-axis (ranging from 105 to 100 eV) and angle in degrees on the y-axis (ranging from 30 to 80 degrees). The plot displays several peaks that shift and change intensity as the angle of detection is varied, providing information about the chemical state and depth of the elements.</p>												

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